

# Locking plate white 1000 ml

Product sheet

If you opt for a 1000 ml MQ dispenser, you want, of course, to enjoy using the dispenser as long as possible. The solid MediQo-line locking plate protects the contents of the dispenser against damage of daily use. If you want to complete your dispenser, you choose the aluminum mediQo-line locking plate. Mounting is simple you put the plate on the front of the dispenser and lock it with the key supplied. The special window allows for easy control of the contents of the dispenser. If you make your 1000 ml MQ dispenser last, choose the aluminum MediQo-line locking plate. Ask for the quote.

[Bekijk het product online](#)



## Technische omschrijving

Locking plate with window for level control. With lock and key (anti-theft of the refills)